

INSTRUMENT FOR IMAGE SCANNING OF IMMUNOBLOT AND WESTERNBLOT STRIPS WITH A DRYING FUNCTION

- Image scanning directly in reaction wells
- Use with different types of trays
- Active drying of strips
- Very good suppression of parasitic reflections
- Competitive price/performance ratio









DYNACAM

Instrument features

- The process of obtaining test results is eased by image scanning the strips directly in the reaction wells of used trays immediately after previous processing.
- It is not necessary to manipulate with the strips between processing and evaluation. It reduces the possibility of making mistakes that are caused by confusing strips during evaluation.
- An universal tray attachment system allows the insertion of different trays into the instrument.
- The instrument contains a system for actively drying the strips using an air stream. This system noticeably reduces the time needed between processing and evaluating these types of strips where the intensity is influenced by the wetness of their surface.
- ◆ The image is scanned strip by strip using a CCD camera.

The construction of the LED lighting, as well as setting of the camera interest area to the bottom of the reaction wells considerably suppresses unwanted reflections which prevent the evaluation software from properly analysing the image.

- Electronic control of camera shutter speed allows reaching optimal brightness of the scanned image.
- The device is controlled with a PC. The software controls the processes of scanning and transmission of image data via functions of DLL library. Communication between PC and the device is carried out via USB interface.
- Universal instrument able to work with different evaluating SWs.
- DYNACAM brings a next step of automation in blot processing when used with a blot processor (for example DYNABLOT Plus).

Technical specifications

Picture taking principle	CCD camera, 5 Mpixel, monochromatic
Resolution	500 DPI
Imaging area	460 x 130 mm
Dominant wave length of LEDs in lighting	525 nm
Strip drying	Ambient temperature air blowing
Connection to PC	2 x USB (instrument, camera)
Power supply	24 V DC (external power supply)
Dimensions	620 mm(W) x 390 mm(D) x 270 mm(H)
Weight	16 kg
Cat. No.	D1014-E



Ver_02_10/2016

PRODUCER

DYNEX TECHNOLOGIES, spol. s r.o.

Place of business: Lidicka 977, 273 43 Bustehrad, The Czech Republic Headquarters: Vodickova 791/41, 110 00 Praha 1, The Czech Republic tel.: + 420 220 303 600, fax: +420 224 320 133 email: office@dynex.cz

www.dynex.cz



DISTRIBUTOR



