



NEXT GENERATION OF AUTOMATED INSTRUMENT USED TO PROCESS THE IMMUNOBLOT AND WESTERNBLOT ASSAYS

MAIN FEATURES

- Walk-away system
- ◆ Automates reagents dispensing, incubating with rocking, strips washing, final drying
 - strips are ready to scan and evaluate
- Capacity up to 50 strips
- ◆ 8 reagent channels fitted with peristaltic pumps
- CE IVD marked

MAJOR IMPROVEMENTS

- Fluid lack sensor in each dispensing channel
- Strip drying function (heated tray and forced airflow)
- Easily accesible and space-saving drawer for reagents
- Improved speed resolution and adjustable angle of tray rocking
- Upgraded incubation timing management
- Low vacuum level alarm
- ◆ Modern colour touch screen for easy and comprehensible operation
- Enhanced assay list capacity to 50 positions
- Great traceability thanks to run logs archiving
- New PC software for comprehensible assay creating (transfer to instrument memory via USB)



OEM version available

OPEN platform

Made in

the Czech

Republic

Complete the workflow: Scan the strips with our DYNACAM instrument.



Ver.01_11/2017



PRODUCER

DYNEX TECHNOLOGIES, spol. s r.o.

Place of business: Lidicka 977, 273 43 Bustehrad, The Czech Republic Headquarters: Vodickova 791/41, 110 00 Praha 1, The Czech Republic tel.: + 420 220 303 600, fax: +420 224 320 133

email: export@dynex.cz www.dynex.eu



